75 Mill Street  
Stoughton MA, 02070 USA.  
TEL.: +1.781.344.3670 | FAX: +1.781.341.9693  
[www.syncmedical.com](http://www.syncmedical.com)

RSQR, LTD.

Ludgate House  
107-111 Fleet Street, London EC4A 2AB  
Tel +44 (0)207-936-9412  
[info@rsqa.co.uk](mailto:info@rsqa.co.uk) [www.rsqa.co.uk](http://www.rsqa.co.uk)

0459

## Ethos Spine™ Cervical Plate System – Package Insert

**Description:** The Ethos Spine™ Cervical Plate System is intended as a temporary internal fixation device used for correction and stabilization of the cervical spine. The system is also intended for aid in the development of a solid spinal fusion.

The Ethos Spine™ System consists of a number of plates and screws of varying lengths. The plates feature lordotic curvature and a transverse plane curvature for an anatomical fit. The set contains ancillary products and instruments. All implant components are machined from implant grade titanium alloy, Ti-6Al-4V (ELI) per ASTM F-136.

### **Indications and Usage:**

The SyncMedical Ethos Spine™ Cervical Plate System is intended for use in the cervical area of the spine. If properly used, the Ethos Spine™ system will aid in the development of solid spinal fusion. The surgeon must be aware of the medical and surgical aspects of the implant, as well as the material limitations of metallic implants. The patient must be instructed in the limitations of the implant, especially with regard to weight bearing and other stresses seen in the body before spinal fusion occurs. The patient must be instructed that non-compliance with postoperative instructions could lead to failure of the device which may require need for additional surgery to remove failed components.

### **Indications for use** are as follows:

The SyncMedical Ethos Spine™ Cervical Plate System is intended for use in anterior cervical fixation for the following indications:

- i. Degenerative Disc Disease (DDD) (defined as neck pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies),
- ii. Spondylolisthesis,
- iii. Trauma (i.e., fracture or dislocation)
- iv. Spinal stenosis,
- v. Deformities or curvatures (i.e., scoliosis, kyphosis, and/or lordosis),
- vi. Tumor,
- vii. Pseudoarthrosis; and
- viii. Previous failed fusion.

The SyncMedical Ethos Spine™ Anterior Cervical Plate System is indicated for stabilizing the cervical spine from C2 to C7.

Unicortical or bicortical fixation from the anterior face of the vertebral bodies may be employed using either the variable angle or fixed angle screws. This product may be employed as an internal fixation device during the time interval required for arthrodesis.

### **Contra-Indications for use** are as follows:

- i. Morbid obesity;
- ii. Mental illness;
- iii. Alcoholism and/or drug abuse;
- iv. Pregnancy;
- v. Metal sensitivity/allergies;
- vi. Severe osteopenia;
- vii. Active infection near the surgical site; and
- viii. Patients unwilling or unable to follow post-operative care instructions.

### **Possible adverse effects:**

- Nonunion.
- Delayed union.
- Bending, fracture, or loosening of implant.
- Metal sensitivity or allergic foreign-body reaction.
- Infection
- Stress shielding
- Pain, discomfort, or other abnormal sensations due to the device.
- Erosion of blood vessels due to the proximity of the device. Such erosion could lead to vessel hemorrhage.
- Bursitis.
- Dural Leak.
- Paralysis
- Death

**WARNING:**

Correct selection of the implant is extremely important. The chance of adequate fixation is increased by the selection of the proper size of the implant. Proper selection of the implant can minimize, but not eliminate risks. No internal fixation device can withstand activity levels equal to those withstood by normal healthy bone. No implant can be expected to withstand unsupported weight-bearing stress indefinitely.

Mixing metals can cause corrosion. Titanium and stainless steel in contact will accelerate the corrosion process in the stainless steel. The presence of corrosion may accelerate the fatigue failure of the implanted devices. The amount of metallic compounds released into the body will be increased. Therefore, internal fixation devices which may come into contact with other metal devices must be made from similar or compatible materials.

Implants may fail when subjected to the increased loading associated with delayed union or nonunion. Internal fixation devices are intended as load-sharing devices until normal healing occurs. If normal healing does not occur or is delayed, the device may break due to fatigue. Notches, scratches, or bending of the implant during the course of surgery may contribute to early failure. Patients should be informed of the risks of implant failure.

**Patient Selection** the following factors may be of importance when selecting patients for internal fixation devices:

- Certain degenerative diseases
- Foreign body sensitivity
- Senility, mental illness, alcoholism, or drug abuse.
- Obesity
- The patient's occupation or activity.

**PRECAUTIONS**

- **Explanted surgical implants must never be reused. This product is single use only.** Even though a device appears undamaged, it may have defects and internal stress patterns which may lead to premature failure.
- **Correct handling of the implants is extremely important.** Surface damage or alterations may produce defects in surface finish and may create internal stresses that could facilitate the eventual breakage of the implant.
- **The patient must be adequately instructed in the appropriate postoperative care** The patient's ability and willingness to follow instructions are some of the most critical aspects of successful bone healing. The patient must be instructed in the limitations of the implant, and that physical activity and full weight bearing may cause premature failure of internal fixation devices by loosening, bending, or fracture. The patient must be instructed that an implant does not have the strength of normal healthy bone, and that device failure may occur if excessive loading is placed on the implant.
- **Care must be taken in removal of the implant after healing.** Particularly in young active patient, implants may loosen, fracture, corrode, migrate, increase the risk of infection, cause pain, or stress shield bone – even after normal healing. The surgeon should consider the risks and benefits when deciding whether or not to remove an implant. Implant removal should be followed by careful postoperative management to avoid re-fracture. For certain patients, the surgeon may elect to not remove the implant in order to eliminate the risks of another surgery.

**POSTOPERATIVE MOBILIZATION**

Careful patient handling for the first two to four months following the operation is very important while the fusion mass matures and becomes able to share load with the implant. Instructing the patient how to reduce stress on the implants is an important part of the effort to avoid clinical problems that may accompany failed fixation.

External immobilization is recommended until the maturation of the fusion mass is confirmed radiographically.

Once the fusion mass has healed, the device may be removed to prevent stress-shielding of the fused bone. As in all patient care issues, this is left to the discretion of the operating surgeon.

**CLEANING AND STERILIZATION**

Implants and instruments are supplied clean but NOT STERILE. They must be sterilized by the hospital prior to use. Remove all packaging materials prior to sterilization. For a  $10^{-6}$  Sterility Assurance Level, the Ethos Spine™ Cervical Plate System implants and instruments are recommended to be steam sterilized by the hospital using one of the two sets of process parameters listed below.

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME
Steam	Pre-Vacuum (Wrapped)	273°F to 275°F (134°C to 135°C)	10 Minutes
Steam	Gravity Displacement (Wrapped)	270°F (132°C)	30 Minutes

**NOTE:** Because the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g. temperatures, times) used for equipment.

**CAUTION:** Federal Law (U.S.A) restricts these devices to sale by or on order of a physician.

75 Mill Street  
Stoughton MA, 02070 Estados Unidos  
TEL. +1.781.344.3670 | FAX: +1.781.341.9693  
[www.syncmedical.com](http://www.syncmedical.com)

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